



Ideal for Sanitizing All Types of Non-Critical Items

Caddy-UV - 1



Caddy-UV - 2

The Caddy-UV – 1 is equipped with a stainless steel storage tray. This sanitizer has one shelf and one powerful UV light on top for good exposure. This unit comes with highly reflective walls.

Outside Dimensions: 14-3/4"W x 8-3/8"D x 9-1/2"H.

Inside Dim.: 13-3/4"W x 7-1/2"D x 7-1/4" H

Volts: 120 volts

Watts: 11 watts

Wt.: 14 lbs.

Price: \$395

Warranty: 1 year



- Cell phones, tablets
- Stethoscopes
- Thermometers
- Blood pressure cuffs
- Oximeters
- Cell phones & Tablets
- MARs scanners
- Penlights
- Scissors
- Otoscopes
- Computer tablets
- Keys/ID cards
- Writing instruments
- Eyeglasses
- All desktop items

The Caddy-UV – 2 is equipped with a timer. This sanitizer has two shelves and two powerful UV lights, one on top and one on bottom for maximum exposure. This unit comes with highly reflective walls.

Outside Dimensions: 15.5" H x 16" D x 16.2" W

Inside Dimensions: 8" H x 9" D x 15" W

Volts: 120 volts

Watts: 20 watts

Wt.: 25 lbs.

Price: \$495

Warranty: 1 year



MRSA-UV™
 Mobile Room Sanitizing Agency
 561-531-9398