



# EMS-UV 24/7 Air Sanitizer





**Advanced UV Office, Home and Ambulance Air Sanitizer**

*Ideal for Hospitals, Clinics, Schools, Ambulances, Day Care, Transport Vehicles, Assisted Living, Restrooms...*

**EMS-UV™** is the most advanced 24/7 air sanitizer on the market. Utilizing ultraviolet lamp technology, this product kills up to 99.9% of airborne bacteria, viruses, black mold and fungal growth, reduces odors and respiratory allergens and reduces airborne transmission of influenza, colds, and viruses.

- Wall mount or table top
- Self-cleaning with no filters to replace
- Rust proof all aluminum alloy profile
- Strong 15w germicidal UV-C lamp
- Nano coated titanium dioxide PCO plates
- Breaks down viruses, air pollutants and VOC
- Negative ions effect on dust removal
- Reflective aluminum chamber
- Quiet operation. Rooms can be fully occupied
- Low power consumption
- Detachable power cord
- 13.8" x 5.1" x 4" (350 x 130 x 100mm)
- Wt. 4.4 lbs. (2 kg.)
- Long lasting and inexpensive UV-C bulb

**EMS-UV™** is ideal for rooms 500 sq. ft. or less and is designed to run 24/7. Bulb lasts 8,000 hours.



**MRSA-UV™**  
Mobile Room Sanitizing Agency  
561-531-9398

