



OBELISK™ UV ROOM SANITIZER



Our Most Affordable, State-of-the-Art, Portable UVC Room Sanitizer

Our advanced **Obelisk™ UV Total Room Sanitizer** is the most affordable, lightweight, portable and extremely powerful UV surface and air sanitizer. Using eight (8) concentrated high-energy UV-C lamps in each unit, this device is capable of disinfecting average size room very quickly and even **faster** when used in tandem.

The Obelisk™ UV Sanitizer can be easily carried from one room to the next while the lamps are always protected by our high quality stainless steel cover.

The Obelisk™ also comes with two ultra-sensitive bacterial verification test kits and our mobile hand truck and stand which is essential when elevated heights are required. Rooms with beds are best treated when the Obelisk™ is used in tandem.


Our Obelisk™ is capable of safely disinfecting environs & air in a room using UVC energy in as little as a 15 minutes. UVC has been long proven in independent laboratory testing to successfully disinfect against all microbes including the most difficult bacteria and viruses. Having eight (8) powerful, shatterproof coated UV-C Germicidal bulbs available allows for quick environmental and air sanitizing.

Your Obelisk™ comes with everything necessary to go to work right away. It includes UV glasses, remote control, Laser Tape Measure, UV Dose Chart, quick instructions, full instructions and our tablet/phone app to calculate exposure times.



Should only be used only as an adjunct to your existing cleaning program.





Model:	40-0633A	40-0673A
Volts:	120v	220-240v
Amps:	4.5A	2.5A
Hertz:	50/60Hz	50/60Hz
Maximum Treatment Time:	24 Hours	
Weight:	31 lbs.	
Dimensions:	Length: (Inches)	13-15/16"
	Width: (Inches)	12-1/16"
	Height: (Inches)	31-9/16"
Lamp Data:	Lamp Model:	05-0892
	Number of Lamps:	8
	Power Consumption: Ⓞ (Watts/Lamp)	68.75W (550W Total)
	Ultraviolet Output: Ⓢ (Watts/Lamp)	17.8W (143W Total)
	Average Effective Life: (Hours)	13,000