



Transport-UV™ is a fixed, direct, surface to air, interior cabin sanitizer. It takes only 15-30 min. to completely disinfect any compartment.



Our products always come with our unique rapid bacteria verification test kit.

TRANSPORT-UV™ RAPIDLY DISINFECTS AIR AND SURFACES IN JUST MINUTES!!

Our rapid **Transport-UV™** is designed to quickly disinfect the air and surfaces within all transport vehicles. It is a fixed unit that is securely mounted to the ceiling of cabin. **Transport-UV™** takes only 10 to 30 minutes to thoroughly sanitize a cab and all of the air and surfaces within it.


The **Transport-UV™** is available with either a safety wall switch or remote control on/off switch. All lamps have been made shatterproof. Additionally, you may choose between 1 to 4 UV lamps in each unit as well as lengths from 18" to 64". The more lamps you have, the faster the treatment time.





Because ambulances and other transport vehicles have a relatively small cabin space, it takes only a few minutes to thoroughly disinfect all surfaces. In addition, UVC also destroys VOC's and other odors.

All forms of transport vehicles have been implicated in numerous infections including MRSA, TB, CRE, C. diff and even Norovirus to name just a few. Both responders and patients are at risk of contracting severe multi-drug resistant diseases.

Recent studies reveal that only 50% of surfaces within a given vehicle are disinfected during cleaning. More importantly, several prior studies show that inadequate cleaning of vehicles occupied by patients with MDR (multidrug resistant) pathogens places the next patient and EMS personnel in these cabs at an increased risk of acquiring these organisms.



ORDERING INFORMATION				
Part No.	Desc.	Plug Type	Price	
12-1L = 12 inches long, 1 lamp		A, B, C, D, E, F, G, H, I		
B133-12-1L 12"	UV LIGHT:L-1 3.1W		\$625.00/each	
B133-12-2L	UV LIGHT:L-2 6.2W		\$655.00/each	
B133-12-4L	UV LIGHT:L-4 12.4W		\$785.00/each	
B133-18-1L 18"	UV LIGHT:L-1 5.8W		\$645.00/each	
B133-18-2L	UV LIGHT:L-2 11.6W		\$685.00/each	
B133-18-4L	UV LIGHT:L-4 23.2W		\$935.00/each	
B133-24-1L 24"	UV LIGHT:L-1 8.5W		\$670.00/each	
B133-24-2L	UV LIGHT:L-2 17W		\$725.00/each	
B133-24-4L	UV LIGHT:L-4 34W		\$975.00/each	
B133-30-1L 30"	UV LIGHT: L-1 30"		\$685.00/each	
B133-30-2L	UV LIGHT: L-2 30"		\$855.00/each	
B133-30-4L	UV LIGHT: L-4 30"		\$1,100/each	
B133-36-1R 36"	UV LIGHT:L-1 13.8W		\$800/each	
B133-36-2L	UV LIGHT:L-2 26W		\$865.00/each	
B133-36-4L	UV LIGHT:L-4 52W		\$1,150/each	
B133-48-1L 48"	UV LIGHT:L-1 19.3W		\$795.00/each	
B133-48-2L	UV LIGHT:L-2 38.6W		\$915.00/each	
B133-48-4L	UV LIGHT:L-4 77.2W		\$1,175.00/each	
B133-64-1L 64"	UV LIGHT:L-1 25W		\$1,100.00/each	
B133-64-2L	UV LIGHT:L-2 50W		\$1,225.00/each	
B133-64-4L	UV LIGHT:L-4 100W		\$1,995.00/each	



The **Transport-UV GERMICIDAL UV FIXTURE** series are stainless steel ceiling or wall mounted fixtures that irradiate room air and exposed surfaces with direct germicidal ultraviolet rays. This treatment is very effective due to the long unobstructed path of the rays. Germicidal Ultraviolet Fixtures should be utilized during periods when the area is unoccupied or when occupants take appropriate protective measures. Germicidal Ultraviolet Fixtures are open fixtures that irradiate room air and exposed surfaces with direct germicidal rays. Optional Safety Wall Switch or Remote Control.

FEATURES:

Transport-UV™ Slim-line base, single pin, lamp holder set (includes one stationary end and one spring end) holds lamp securely while providing convenient, "quick and easy" lamp change. Includes protective grating. Order with 1 to 4 lamps.

Transport-UV™ State-of-the-art electronic ballast specifically developed for the operation of ultraviolet lamps. These versatile, instant start ballasts provide higher lamp output; are smaller, lighter, more efficient and operate cooler than magnetic ballasts.

Standard items are 120v 50/60Hz with shatterproof UV bulbs and protective mesh. For international customers, a 220v 50/60Hz item may be special ordered for a slightly higher price; please call us for a quote. Optional 12 volt DC inverter.

