



Our products always come with our sample pack of UV-C verification test labels.

TRANSPORT-UV™ RAPIDLY DISINFECTS AIR AND SURFACES IN JUST MINUTES!!

Our Slim-Lined **Transport-UV™** is designed to quickly disinfect the air and surfaces within all transport vehicles. It is a fixed unit that is securely mounted to the ceiling of cabin. **Transport-UV™** takes only 10 to 30 minutes to thoroughly sanitize a cab and all of the air and surfaces within it.

The **Transport-UV™** is comes with our advanced multi-featured remote control on/off switch so that the unit can be easily turned on from outside of the cab as well as setting the treatment time.

The **Transport-UV** unit comes with twin types of motion-sensor shut-off systems, PIR and radar. This insures enormous safety should someone inadvertently enter the vehicle. The Transport-UV unit has twin UV-C germicidal lamps and are three feet in length.

Because ambulances and other transport vehicles have a relatively small cabin space, it takes only a few minutes to thoroughly disinfect all surfaces. In addition, UVC also destroys VOC's and other odors.

All forms of transport vehicles have been implicated in numerous infections including MRSA, TB, CRE, C. diff and even Norovirus to name just a few. Both responders and patients are at risk of contracting severe multi-drug resistant diseases.

Recent studies reveal that only 50% of surfaces within a given vehicle are disinfected during cleaning. More importantly, several prior studies show that inadequate cleaning of vehicles occupied by patients with MDR (multidrug resistant) pathogens places the next patient and EMS personnel in these cabs at an increased risk of acquiring these organisms.

